

# MINNCARE®

## Peracetic Acid Test Strips



Minncare 1% Test Strips



### Minncare Test Strips Offer:

- Fast test results
- Convenient ready-to-use test system
- No mixing or dilutions
- Simple, easy to use procedures



Minncare Residual Test Strips

# Quick results and easy-to-read indicators.

Minnicare® test strips provide quick results with easy-to-read indicators. Use the Minnicare® 1% TS test strips as a pass/fail measurement for adequate concentration of sterilant after dilution. Then, verify residual levels after rinse-out using Minnicare® Residual Test Strips.

Minnicare Peracetic Acid Test Strips represent the leading edge of test strip technology. Designed by the producers of Minnicare® Cold Sterilant, these test strips combine ease-of-use with precise results.

## Minnicare® 1% Test Strips

Minnicare® 1% Test Strips verify adequate concentration of a 1% Minnicare® solution during sanitization of your high purity water system. These convenient dip-and-read strips can be used at any test port in your system to provide an accurate pass/fail measurement of Minnicare® solution.

## Minnicare® Residual Test Strips

Minnicare® Residual Test Strips are easy-to-use and are conveniently packaged 100 test strips per container. Each strip has a "reaction zone" pad located at one end of the strip. When stored and used properly, the test strip pad can indicate residual Minnicare® levels as low as 1 part per million (ppm). The reaction zones are set to provide six different reading levels at 100, 30, 10, 3, 1 and 0 (zero) parts per million.

Both the Minnicare® 1% TS and the Minnicare Residual Test Strips may be used up to the date of expiration printed on the container even after opening. There is no need to separately record and verify container opening dates.

**To Order:**

**Call: (800) 328-3370**

### Ordering Information

Model	Description
78338-000 . . . .	Minnicare® Residual Test Strips (100 Strips)
89330-000 . . . .	Minnicare® 1% TS Peracetic Acid Test Strips (100 Strips)



Minntech Corporation  
14605 28th Avenue North  
Minneapolis, MN 55447-4822  
U.S.A.  
Tel: (763) 553-3300  
Toll Free: (800) 328-3370  
Fax: (763) 553-3387

Minntech B.V.  
Sourethweg 11  
6422PC Heerlen  
The Netherlands  
Tel: (31) 45 5 471471  
Fax: (31) 45 5 429695  
www.minntech.com

Minntech Japan Corporation  
4F Sumitomo Seimei Yochomachi Bldg.  
10-10 Yochomachi, Shinjuku-ku  
Tokyo 162-0055  
Japan  
Tel: (81) 3 (3225) 8680  
Fax: (81) 3 (3225) 8681